



CU Index™
Functional Anti-FcεR
CU Autoantibody Signature™
Auto Immune Chronic Urticaria Test
Test Code: 2103

CPT CODE	86352
SPECIMEN REQUIREMENT	
Type:	Serum
Volume:	1 mL
Handling:	Blood should be collected and allowed to clot prior to centrifugation: (NOTE: Patients taking calcineurin inhibitors should stop their medication for 72 hours prior to draw. Patients taking prednisone should be off their medication for 2 weeks prior to draw.
Storage:	Samples can be stored for 1 week refrigerated, if the specimen is to be held for more than a week, it should be stored frozen until shipped to the lab.
Shipping:	The specimen can be shipped via overnight courier at ambient temperature. If the specimen is to be held for more than a week it should be stored frozen until shipped to the lab. NOTE: Serum is the ONLY acceptable specimen.
METHODOLOGY	
SET-UP SCHEDULE	Monday thru Friday
STANDARD TURNAROUND TIME	2-3 days from receipt of specimen