



## BasoFunction™ HRT

Test Code: See Listing Below

<b>CPT CODE</b>	Basophil Activation Test: 86343, 83088 Each Allergen: 86343, 83088
<b>SPECIMEN REQUIREMENT</b>	Type: Whole Blood Volume: 8-10 mL (for 5 antigens) Handling: Sodium heparin tube. Keep blood at room temperature at all times. <u>Do Not Use Lithium Heparin tube.</u> Storage: Ambient Shipping: Ship <b>OVERNIGHT</b> . Samples are received <b>Monday thru Friday ONLY</b> . Call lab with tracking number. Specimen <b>MUST BE</b> received within 48 hours of draw.
<b>METHODOLOGY</b>	Cell Stimulation/ELISA
<b>SET-UP SCHEDULE</b>	Monday through Friday
<b>STANDARD TURNAROUND TIME</b>	3-5 business days from receipt of specimen
<b>CLINICAL UTILITY</b>	The BasoFunction test may be useful with patients who have a suggestive history and for whom skin testing is not appropriate or available.
<b>STABILITY</b>	Not stable

Analgesic		Anesthetic		Antibiotic		Food Additives	
501361	Acetaminophen	501375	Atracurium	501381	Amoxicillin	501366	MSG
501369	Diclofenac	501353	Lidocaine	501382	Ampicillin	501359	Potassium Metabisulphite
501370	Dipyron/ Metamizole	501376	Mivacurium	501383	Cefamandole	501364	Tartrazine
501360	Ibuprofen	501377	Pancuronium	501384	Cefazolin	501362	Quinoline Yellow
501371	Indomethacin	501367	Propofol	501385	Cefuroxime	501363	Indigo Carmine
501365	Lys-Aspirin	501378	Rocuronium	501351	Cephalosporin C	501379	Sodium Benzoate
501372	Mefanamic Acid	501357	Suxamethonium	501352	Ciprofloxacin	501394	Sodium Nitrite
501373	Phenylbutazone Acid	501380	Vecuronium	501355	Penicillin: major	501368	Sodium Salicylate
501374	Propyphenazone			501354	Penicillin: minor	501387	Sunset Yellow FCF
				501356	Sulfamethoxazole	501389	Chromotrope B
				501358	Tetracycline	501390	Amaranth
				501386	Trimethoprim	501395	New Coccine
						501391	Erythrosine
						501388	Patent Blue V
						501392	Brilliant Black BN